

PATHOPAK



**SAMPLE
TRANSPORT
PACKAGING**



PathoPak sample transport packaging solutions are strong, durable, reusable, rigid secondary containers compliant with IATA, ADR, and CFR 49 (DOT) transport regulations for UN Class 6.2 Infectious Substances Category A and UN3373 Category B Infectious Substances.

PathoPak packaging is suitable for shipping a wide range of primaries, including blood tubes, specimen containers, swabs, blood collection/transfusion bags and medical devices.

These reusable rigid systems are cost-effective and are available in four sizes: 800mL, 1L, 2L and 3L, allowing you to choose the right size for your payload, minimising wasted space.

THIS PRODUCT INCLUDES:

- 95kPa Tested Secondary Container
- Absorbent
- Bubble wrap
- Labels and UN markings

THIS PRODUCT CAN BE USED FOR:

- Infectious samples
- Diagnostic samples
- Clinical trials
- Screening programmes
- Eradication schemes

FEATURES:

- Reusable
- Ergonomically engineered bottle design for safe manual handling
- 4GU certified
- Suitable for cleaning by autoclave
- Full refurbishment services available

For more information please contact Intelsius UK cs@intelsius.com



PATHOPAK RANGE

Product Code	Product Name	Case Quantity	External Dimensions	Payload Dimensions (Bottle opening diameter x internal height)	Payload Volume	Mass	Volumetric Weight L/6
PP034	DGP PathoPak 800mL with 6 Bay Absorbent	20	112 x 110 x 150 mm 4.4 x 4.3 x 5.9 in	90 x 134 mm 3.5 x 5.3 in	800mL	0.19 kg 0.41 lbs	0.31 kgvol ^c 0.68 lbsvol ^c
PP011	DGP PathoPak 1L with 50ml Absorbent	20	112 x 110 x 190 mm 4.4 x 4.3 x 7.5 in	90 x 170 mm 3.5 x 6.7 in	1 L	0.24 kg 0.53 lbs	0.39 kgvol ^c 0.86 lbsvol ^c
PP012	DGP PathoPak 2L with 6 Bay Absorbent	10	153 x 156 x 193 mm 6.0 x 6.1 x 7.6 in	125 x 170 mm 4.9 x 6.7 in	2 L	0.45 kg 0.99 lbs	0.77 kgvol ^c 1.69 lbsvol ^c
PP013	DGP PathoPak 3L with 6 Bay Absorbent	10	160 x 160 x 280 mm 6.3 x 6.3 x 11.0 in	125 x 255 mm 4.9 x 10.0 in	3 L	0.6 kg 1.32 lbs	1.19 kgvol ^c 2.63 lbsvol ^c

UNITED KINGDOM: DGP Intelsius Ltd
1 Harrier Court, Airfield Business Park,
Elvington, York, YO41 4EA,
United Kingdom

t: +44 (0) 1904 607 390
f: +44 (0) 1904 607 391
cs@intelsius.com

AMERICA: DGP Intelsius LLC
1399 Perry Rd, Suite 190,
Plainfield, IN 46168,
United States of America

t: +1 317 452 4006
f: +1 317 290 4621
americasorders@intelsius.com

MALAYSIA: DGP Intelsius Sdn Bhd
Plot 11, Lorong Bemban 2, Estet
Perindustrian Bemban, 31000
Batu Gajah, Malaysia

t: +605 3666 399
f: +605 3669 399
compliance@intelsius.com

GERMANY: DGP Intelsius GmbH
Eulerweg 11
64291 Darmstadt
Germany

t: +49 6151 4936650
f: +49 (0) 6151 4936651
info@intelsius.de

IRELAND: DGP Intelsius Ireland Ltd.
Unit 8, Coolcaslagh Industrial Estate,
Coolcaslagh, Killarney
Co. Kerry, V93 K006, Ireland

t: +353 (0)64 668 61 71
salesireland@intelsius.com

INDIA: DGP Intelsius Cold Chain Pvt. Ltd.
Sun Srishti, C 10, Saki Vihar Road,
Chandivali Powai. Mumbai - 400072
Maharashtra - India

+91 22 49746536
compliance@intelsius.com